

Product Datasheet

CD289 antibody (orb430568)

Description	Goat polyclonal antibody to CD289
Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic peptide LSLKYNNLTQVPRQLPPSLEY-C corresponding to amino acids 204-224 within the N-terminal region of mouse CD289.
Target	CD289
Preservatives	0.1% Sodium Azide (NaN ₃)0.1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	2.0mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Application notes	(N-TERMINAL)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.