

Product Datasheet

ATG5-ATG12 COMPLEX antibody (orb430570)

Description	Rabbit polyclonal antibody to ATG5-ATG12 COMPLEX
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide corresponding to amino acids 262-275 of mouse Atg5. The corresponding sequence differs by one amino acid in human.
Target	ATG5-ATG12 COMPLEX
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline.
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(C-TERMINAL)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.