

Product Datasheet

DELTA-LIKE PROTEIN 3 antibody [Out of stock] (orb430572)

Description	Rabbit polyclonal antibody to DELTA-LIKE PROTEIN 3
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	DELTA-LIKE PROTEIN 3
Concentration	1.0mg/ml
Note	For research use only
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.