

Product Datasheet

IL-1 ALPHA antibody (orb430584)

Description	Goat polyclonal antibody to IL-1 ALPHA
Species/Host	Goat
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FA, WB
Immunogen	Recombinant rat IL-1 alpha.
Target	IL-1 ALPHA
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - lyophilized
Concentration	1.0 mg/ml
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.