

Product Datasheet

VEGF antibody (Biotin) (orb430600)

Description	Rabbit polyclonal antibody against VEGF conjugated to Biotin
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Recombinant rat VEGF.
Target	VEGF
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - lyophilized
Concentration	0.1mg/ml
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.