

Product Datasheet

STEM CELL FACTOR antibody (orb430602)

Description	Rabbit polyclonal antibody to STEM CELL FACTOR
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FA, WB
Immunogen	Recombinant rat SCF.
Target	STEM CELL FACTOR
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - lyophilized
Concentration	1.0mg/ml
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.