

Product Datasheet

C1q antibody (FITC) (orb430606)

Description	Sheep polyclonal antibody against C1q conjugated to FITC
Species/Host	Sheep
Reactivity	Human
Conjugation	FITC
Tested Applications	IH-Fr
Immunogen	Human C1q, purified from plasma.
Target	C1q
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	10 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.