

## Product Datasheet

### VON WILLEBRAND FACTOR antibody (orb430610)

|                            |                                                                                                   |
|----------------------------|---------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Sheep polyclonal antibody to VON WILLEBRAND FACTOR                                                |
| <b>Species/Host</b>        | Sheep                                                                                             |
| <b>Reactivity</b>          | Human                                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                                      |
| <b>Tested Applications</b> | ELISA, FACS, ICC, IF, IHC-Fr                                                                      |
| <b>Immunogen</b>           | Purified human von Willebrand factor prepared from citrated human plasma                          |
| <b>Target</b>              | VON WILLEBRAND FACTOR                                                                             |
| <b>Preservatives</b>       | 0.1% Sodium Azide (NaN <sub>3</sub> )0.1% EACA0.01% Benzamidine1 mM EDTA. Glycine buffered saline |
| <b>Form/Appearance</b>     | Purified IgG - liquid                                                                             |
| <b>Storage</b>             | +4°C, -20°C if preferred                                                                          |
| <b>Note</b>                | For research use only                                                                             |
| <b>Isotype</b>             | Polyclonal IgG                                                                                    |
| <b>Clonality</b>           | Polyclonal                                                                                        |
| <b>Purity</b>              | Purified                                                                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                   |