

Product Datasheet

TSC1 antibody (orb430614)

Description	Rabbit polyclonal antibody to TSC1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IH-Fr, IHC-P, WB
Immunogen	A 15 amino acid peptide located near the centre of human TSC1.
Target	TSC1
Preservatives	0.02% Sodium Azide. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.