

## Product Datasheet

### ADIPONECTIN antibody (orb430617)

|                            |                                                                          |
|----------------------------|--------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to ADIPONECTIN                                |
| <b>Species/Host</b>        | Rabbit                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                             |
| <b>Tested Applications</b> | FC, IHC-P, WB                                                            |
| <b>Immunogen</b>           | A 15 amino acid peptide located near human Adiponectin carboxy-terminus. |
| <b>Target</b>              | ADIPONECTIN                                                              |
| <b>Preservatives</b>       | 0.02% Sodium Azide (NaN <sub>3</sub> ). Phosphate buffered saline        |
| <b>Form/Appearance</b>     | Purified IgG - liquid                                                    |
| <b>Concentration</b>       | 1.0mg/ml                                                                 |
| <b>Storage</b>             | +4°C, -20°C if preferred                                                 |
| <b>Note</b>                | For research use only                                                    |
| <b>Application notes</b>   | (C-TERMINAL)                                                             |
| <b>Isotype</b>             | Polyclonal IgG                                                           |
| <b>Clonality</b>           | Polyclonal                                                               |
| <b>Purity</b>              | Purified                                                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                          |