

Product Datasheet

FGF BASIC antibody (orb430646)

Description	Rabbit polyclonal antibody to FGF BASIC
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FA, IHC-P, WB
Target	FGF BASIC
Preservatives	None present.. Phosphate buffered saline
Form/Appearance	Purified IgG - lyophilized
Concentration	1.0 mg/ml
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.