

Product Datasheet

SYK antibody (orb430669)

Description	Goat polyclonal antibody to SYK
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide with sequence C-KYLQQNRHVKDKN from the internal region of the SYK sequence.
Target	SYK
Preservatives	0.02% Sodium Azide (NaN ₃) 0.5% Bovine Serum Albumin. TRIS buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.