

## Product Datasheet

### DYKDDDDK TAG antibody (orb430674)

|                            |                                                                                                                                                                             |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to DYKDDDDK TAG                                                                                                                                  |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                      |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                |
| <b>Tested Applications</b> | ELISA, WB                                                                                                                                                                   |
| <b>Immunogen</b>           | Keyhole limpet haemocyanin conjugated Enterokinase cleavage site peptide DYKDDDDK. Residues of cysteine and glycine were used to facilitate coupling at the C-terminal end. |
| <b>Target</b>              | DYKDDDDK TAG                                                                                                                                                                |
| <b>Preservatives</b>       | 0.1% Sodium Azide (NaN <sub>3</sub> ). Phosphate buffered saline                                                                                                            |
| <b>Form/Appearance</b>     | Purified IgG - liquid                                                                                                                                                       |
| <b>Concentration</b>       | 0.5 mg/ml                                                                                                                                                                   |
| <b>Storage</b>             | +4°C, -20°C if preferred                                                                                                                                                    |
| <b>Note</b>                | For research use only                                                                                                                                                       |
| <b>Isotype</b>             | Polyclonal IgG                                                                                                                                                              |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                  |
| <b>Purity</b>              | Purified                                                                                                                                                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                             |