

Product Datasheet

SNAP-25 antibody (orb430701)

Description	Goat polyclonal antibody to SNAP-25
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Peptide sequence C-DEANQRATKMLGSG corresponding to the C-terminus of SNAP-25 (NP_003072.2, NP_570824.1).
Target	SNAP-25
Preservatives	0.02% Sodium Azide 0.5% Bovine Serum Albumin. TRIS buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(C-TERMINAL)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.