

## Product Datasheet

### TIM-1 antibody (orb430744)

|                            |                                                                |
|----------------------------|----------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to TIM-1                            |
| <b>Species/Host</b>        | Rabbit                                                         |
| <b>Reactivity</b>          | Human                                                          |
| <b>Conjugation</b>         | Unconjugated                                                   |
| <b>Tested Applications</b> | WB                                                             |
| <b>Immunogen</b>           | 16 amino acid peptide sequence near the centre of human TIM-1. |
| <b>Target</b>              | TIM-1                                                          |
| <b>Preservatives</b>       | 0.02% Sodium Azide. Phosphate buffered saline                  |
| <b>Form/Appearance</b>     | Purified IgG - liquid                                          |
| <b>Concentration</b>       | 1.0mg/ml                                                       |
| <b>Storage</b>             | +4°C, -20°C if preferred                                       |
| <b>Note</b>                | For research use only                                          |
| <b>Isotype</b>             | Polyclonal IgG                                                 |
| <b>Clonality</b>           | Polyclonal                                                     |
| <b>Purity</b>              | Purified                                                       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                |