

Product Datasheet

ASK1 (phospho-Ser83) antibody (orb430781)

Description	Rabbit polyclonal antibody to ASK1 (phospho-Ser83)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated peptide corresponding to amino acids 76-87 of human ASK-1 protein.
Target	ASK1 (pSer83)
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	2.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.