

Product Datasheet

AURORA-B KINASE antibody (orb430794)

Description	Rabbit polyclonal antibody to AURORA-B KINASE
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide corresponding to amino acids 227-238 of human Aurora-B kinase.
Target	AURORA-B KINASE
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.2mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.