

Product Datasheet

PPAR GAMMA antibody (orb430797)

Description	Rabbit polyclonal antibody to PPAR GAMMA
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Synthetic peptide corresponding to amino acids 255 - 268 of human PPAR gamma isoform 1.
Target	PPAR GAMMA
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.