

Product Datasheet

PPAR ALPHA antibody (orb430799)

Description	Rabbit polyclonal antibody to PPAR ALPHA
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide corresponding to amino acids 1-18 of human PPAR alpha.
Target	PPAR ALPHA
Preservatives	0.1% Sodium Azide (NaN ₃). 0.02 M Potassium Phosphate0.15 M Sodium Chloride
Form/Appearance	Purified IgG - liquid
Concentration	0.98mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(N-TERMINAL)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.