

Product Datasheet

PPAR DELTA antibody (orb430800)

Description	Rabbit polyclonal antibody to PPAR DELTA
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Syntetic peptide corresponding to amino acids 1-14 of mouse PPAR delta.
Target	PPAR DELTA
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(N-TERMINAL)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.