

Product Datasheet

MEK 5 (phospho-Ser311/pThr315) antibody (orb430829)

Description	Rabbit polyclonal antibody to MEK 5 (phospho-Ser311/pThr315)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide corresponding to an amino acid sequence within MEK5 , which includes the phosphorylated residues Ser311 and Thr315.
Target	MEK 5 (pSer311/pThr315)
Preservatives	0.09% Sodium Azide (NaN ₃)0.01% Bovine Serum Albumin50% Glycerol. 10mM HEPES pH 7.5, 150mM NaCl
Form/Appearance	Purified IgG - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.