

Product Datasheet

GLYCOGEN PHOSPHORYLASE BB antibody (orb430838)

Description	Rabbit polyclonal antibody to GLYCOGEN PHOSPHORYLASE BB
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic peptide: internal domain of human brain glycogen phosphorylase, conjugated to Keyhole Limpet Hemocyanin.
Target	GLYCOGEN PHOSPHORYLASE BB
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.