

Product Datasheet

MUSCARINIC RECEPTOR M3 antibody (orb430859)

Description	Rabbit polyclonal antibody to MUSCARINIC RECEPTOR M3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic peptide corresponding to the 3rd cytoplasmic domain of human muscarinic receptor M3, conjugated to Keyhole Limpet Hemocyanin (KLH)
Target	MUSCARINIC RECEPTOR M3
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.