

Product Datasheet

CASPASE-2 antibody (orb430885)

Description	Rabbit polyclonal antibody to CASPASE-2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic peptide corresponding to 6 amino acid sequence from the amino-terminus of human Caspase-2.
Target	CASPASE-2
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.