

Product Datasheet

GALECTIN-3 antibody (Biotin) (orb430907)

Description	Rabbit polyclonal antibody against GALECTIN-3 conjugated to Biotin
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, WB
Target	GALECTIN-3
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - lyophilized
Concentration	0.1mg/ml
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.