

Product Datasheet

ALPHA SYNUCLEIN (phospho-Ser129) antibody (orb430920)

Description	Rabbit polyclonal antibody to ALPHA SYNUCLEIN (phospho-Ser129)
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide corresponding to amino acid sequence surrounding phosphorylated serine 129 of Alpha Synuclein.
Target	ALPHA SYNUCLEIN (pSer129)
Preservatives	0.09% Sodium Azide (NaN ₃) 50% Glycerol 0.01% Bovine Serum Albumin. 10mM Hepes pH 7.5, 150mM NaCl
Form/Appearance	Purified IgG - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.