

Product Datasheet

IL-27 RECEPTOR ALPHA antibody (orb430961)

Description	Rabbit polyclonal antibody to IL-27 RECEPTOR ALPHA
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	A peptide corresponding to amino acids near the carboxy terminus of human IL-27R α precursor.
Target	IL-27 RECEPTOR ALPHA
Preservatives	0.02% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(C-TERMINAL)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.