

Product Datasheet

GAPDH antibody (orb430967)

Description	Rabbit polyclonal antibody to GAPDH
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	A peptide corresponding to a 16 amino acid sequence from near the amino-terminus of GAPDH.
Target	GAPDH
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline.
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(N-TERMINAL)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.