

Product Datasheet

ATG7 antibody (orb430974)

Description	Rabbit polyclonal antibody to ATG7
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IH-Fr, WB
Immunogen	A peptide corresponding to an 17 amino acid sequence from near the carboxy terminus of human ATG7.
Target	ATG7
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(C-Terminal)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.