

Product Datasheet

Bcl-xL antibody (orb431003)

Description	Rabbit polyclonal antibody to Bcl-xL
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A peptide corresponding to 14 amino acids near the carboxy terminus of human Bcl-xL precursor.
Target	Bcl-xL
Preservatives	0.02% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.