

Product Datasheet

DEFENSIN ALPHA 1 antibody (orb431035)

Description	Goat polyclonal antibody to DEFENSIN ALPHA 1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FA, WB
Target	DEFENSIN ALPHA 1
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - lyophilized
Concentration	1.0mg/ml
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.