

Product Datasheet

DEFENSIN BETA 1 antibody (orb431049)

Description	Rabbit polyclonal antibody to DEFENSIN BETA 1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Recombinant human defensin β -1.
Target	DEFENSIN BETA 1
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - lyophilized
Concentration	1.0mg/ml
Storage	-20°C only
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.