

## Product Datasheet

### DEFENSIN BETA 3 antibody (orb431051)

|                            |                                               |
|----------------------------|-----------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to DEFENSIN BETA 3 |
| <b>Species/Host</b>        | Rabbit                                        |
| <b>Reactivity</b>          | Human                                         |
| <b>Conjugation</b>         | Unconjugated                                  |
| <b>Tested Applications</b> | ELISA, WB                                     |
| <b>Target</b>              | DEFENSIN BETA 3                               |
| <b>Preservatives</b>       | None present. Phosphate buffered saline       |
| <b>Form/Appearance</b>     | Purified IgG - lyophilized                    |
| <b>Concentration</b>       | 1.0mg/ml                                      |
| <b>Storage</b>             | -20°C only                                    |
| <b>Note</b>                | For research use only                         |
| <b>Isotype</b>             | Polyclonal IgG                                |
| <b>Clonality</b>           | Polyclonal                                    |
| <b>Purity</b>              | Purified                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.               |