

Product Datasheet

IL-21 RECEPTOR antibody (orb431064)

Description	Goat polyclonal antibody to IL-21 RECEPTOR
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	Synthetic peptide CILEMWNLHPSTLTLTWQDQYEELKDEATS corresponding to amino acids 35-65 within the N-terminal region of human IL-21R.
Target	IL-21 RECEPTOR
Preservatives	0.1% Sodium Azide (NaN ₃) 0.1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Application notes	(N-TERMINAL)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.