

Product Datasheet

ERYTHROPOIETIN RECEPTOR antibody (orb431070)

Description	Sheep polyclonal antibody to ERYTHROPOIETIN RECEPTOR
Species/Host	Sheep
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	Partial cDNA of the extracellular domain of human EPOR expressed as a GST fusion protein in E. Coli.
Target	ERYTHROPOIETIN RECEPTOR
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.