

Product Datasheet

LAT1 antibody (orb431072)

Description	Rabbit polyclonal antibody to LAT1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Keyhole limpet haemocyanin (KLH) conjugated synthetic peptide C-RHRKPEL, representing an internal portion of human LAT1 located towards the C-terminus (aa. 419-425).
Target	LAT1
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Preservative Free
Expiration Date	12 months from date of receipt.