

Product Datasheet

METABOTROPIC GLUTAMATE RECEPTOR 1a antibody (orb431091)

Description	Rabbit polyclonal antibody to METABOTROPIC GLUTAMATE RECEPTOR 1a
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IH-Fr, WB
Immunogen	Peptide corresponding to amino acids 1115-1130 from the C-terminal region of rat mGluR1a.
Target	METABOTROPIC GLUTAMATE RECEPTOR 1a
Preservatives	0.09% Sodium Azide (NaN ₃)0.01% Bovine Serum Albumin50% Glycerol. 10mM HEPES pH 7.5, 150mM NaCl
Form/Appearance	Purified IgG - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.