

Product Datasheet

ATG16 antibody (orb431105)

Description	Rabbit polyclonal antibody to ATG16
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IH-Fr, WB
Immunogen	Synthetic peptide corresponding to an 18 aa sequence from near the centre of human ATG16.
Target	ATG16
Preservatives	0.02% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(INTERNAL REGION)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.