

Product Datasheet

APOLIPOPROTEIN L1 antibody (orb431106)

Description	Goat polyclonal antibody to APOLIPOPROTEIN L1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide C-NNNYKILQADQE from within the C-Terminal region of human APOL1 (NP_003652.2; NP_663318.1).
Target	APOLIPOPROTEIN L1
Preservatives	0.02% Sodium Azide (NaN ₃) 0.5% Bovine Serum Albumin. TRIS buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(C-TERMINAL)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.