

Product Datasheet

CLUSTERIN antibody (orb431167)

Description	Goat polyclonal antibody to CLUSTERIN
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide sequence C-EKALQEYRKKHREE from the C-Terminal region of clusterin (NP_001822.2; NP_976084.1).
Target	CLUSTERIN
Preservatives	0.02% Sodium Azide (NaN ₃) 0.5% Bovine Serum Albumin. TRIS buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(C-TERMINAL)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.