

Product Datasheet

Mcl-1 antibody (orb431169)

Description	Goat polyclonal antibody to Mcl-1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide sequence C-AGVGAGLAYLIR from the C-Terminal region of Mcl-1 (NP_068779.1).
Target	Mcl-1
Preservatives	0.02% Sodium Azide (NaN ₃) 0.5% Bovine Serum Albumin. TRIS buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.5mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(C-TERMINAL)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.