

Product Datasheet

HUMAN IMMUNODEFICIENCY VIRUS 1 gp41 antibody (orb431207)

Description	Goat polyclonal antibody to HUMAN IMMUNODEFICIENCY VIRUS 1 gp41
Species/Host	Goat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthetic peptide CH3CO-YTSLIHSLEESQNQQEKNEQELLELDKWASLWNWF-NH2
Target	HUMAN IMMUNODEFICIENCY VIRUS 1 gp41
Preservatives	0.1% Sodium Azide (NaN3). Phosphate buffered saline
Form/Appearance	Phosphate buffered saline, 0.1% Sodium Azide (NaN3)
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	> 98%
Dilution Range	ELISA: 1:500-10,000, WB: 1:500-10,000
Expiration Date	12 months from date of receipt.