

Product Datasheet

HUMAN IMMUNODEFICIENCY VIRUS 1 gp41 antibody (orb431209)

Description	Goat polyclonal antibody to HUMAN IMMUNODEFICIENCY VIRUS 1 gp41
Species/Host	Goat
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant HIV1 gp41 antigen.
Target	HUMAN IMMUNODEFICIENCY VIRUS 1 gp41
Preservatives	0.1% Sodium Azide (NaN ₃). Isotonic saline solution, pH 6.0-7.0
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.