

Product Datasheet

PAK1 (phospho-Thr84) antibody (orb431239)

Description	Rabbit polyclonal antibody to PAK1 (phospho-Thr84)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Phosphopeptide corresponding to amino acid residues surrounding the phospho-Thr84 (pThr84) region of human PAK1
Target	PAK1 (pThr84)
Preservatives	0.09% Sodium Azide (NaN ₃)1% Bovine Serum Albumin50% Glycerol. HEPES buffered saline
Form/Appearance	Purified IgG - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.