

Product Datasheet

ATF2 (phospho-Ser490/pSer498) antibody (orb431242)

Description	Rabbit polyclonal antibody to ATF2 (phospho-Ser490/pSer498)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Phosphopeptide corresponding to amino acid residues surrounding the phospho-Ser490/Ser498 region of human ATF2
Target	ATF2 (pSer490/pSer498)
Preservatives	0.09% Sodium Azide (NaN ₃)1% Bovine Serum Albumin50% Glycerol. HEPES buffered saline
Form/Appearance	Purified IgG - liquid
Storage	+4°C, -20°C if preferred
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.