

Product Datasheet

Bcl-6 antibody (orb431247)

Description	Rabbit polyclonal antibody to Bcl-6
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, IP, WB
Immunogen	Synthetic peptide mapping to amino acid residues near the C-terminus of human Bcl-6
Target	Bcl-6
Preservatives	30% Glycerol0.5% Bovine Serum Albumin0.01% Thiomersal. Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	0.2 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.