

Product Datasheet

CATALASE antibody (orb431256)

Description	Rabbit polyclonal antibody to CATALASE
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide surrounding amino acid 514 of human catalase
Target	CATALASE
Preservatives	30% Glycerol0.5% Bovine Serum Albumin0.5mM EDTA0.03% Proclin 300. Phosphate buffered saline.
Form/Appearance	Purified IgG - liquid.
Concentration	0.5 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.