

Product Datasheet

CYTOCHROME C antibody (orb431267)

Description	Rabbit polyclonal antibody to CYTOCHROME C
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Synthetic peptide corresponding to residues surrounding amino acid 70 of rat cytochrome c
Target	CYTOCHROME C
Preservatives	30% Glycerol0.5% Bovine Serum Albumin5mM EDTA0.03% Proclin 300. Phosphate buffered saline.
Form/Appearance	Purified IgG - liquid
Concentration	0.2 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.