

Product Datasheet

IP-10 antibody (Biotin) (orb431331)

Description	Rabbit polyclonal antibody against IP-10 conjugated to Biotin
Species/Host	Rabbit
Reactivity	Bovine
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant bovine IP-10 (CXCL10)
Target	IP-10
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Biotin - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.