

Product Datasheet

IL-1 ALPHA antibody (orb431340)

Description	Rabbit polyclonal antibody to IL-1 ALPHA
Species/Host	Rabbit
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant bovine interleukin-1 alpha
Target	IL-1 ALPHA
Preservatives	0.09% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified IgG - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Purity	Purified
Expiration Date	12 months from date of receipt.